

Elite Imaging Table



*Unobstructed 3-D
c-arm imaging*

500lbs/227kg weight capacity

High Quality Hospital Grade Mobile C-Arm Imaging Table

The "Elite" series of imaging tables are designed for high patient throughput utilizing the highest quality components and designed for ease-of-use by the medical staff with patient comfort in mind.

The tabletop has been designed to give maximum unobstructed imaging allowing 360° rotation radiation-free coverage around the longitudinal axis of the carbon fiber tabletop which makes it ideal for 3D imaging studies.

The tabletop pad is also designed to give the patient maximum comfort being manufactured utilizing our patient form-fitting Visco/Memory "Sofcare" High Density foam.

150mm (58") Full 360° Unencumbered Rotational Imaging Area



Model: 0047-TM3

Multi-purpose, 3-movement
Imaging Table

- Motorized Elevation
- $\pm 15^\circ$ Trendelenburg Tilt
- $\pm 15^\circ$ Lateral Tilt

Model: 0047-TM5

Extend-Retract, 5-movement
Imaging Table

- Motorized Elevation
- 15.75" (400mm) Motorized Longitudinal Travel
- $\pm 4"$ (100mm) Motorized Lateral Travel
- $\pm 15^\circ$ Trendelenburg Tilt
- $\pm 15^\circ$ Lateral Tilt





Movement One

- Motorized Elevating Lift
- Lift Travel: (11.75") 300mm
- Optional Travel: (15.75") 400mm



Movement Two

- Motorized Trendelenburg Tilt
- ± 15 Degree Tilt

Movement Three

- Motorized Extend and Retract
- Longitudinal Travel: (15.75") 400mm



Movement Four

- Motorized Lateral Tilt
- Lateral Tilt: ± 15 Degree Tilt



Movement Five

- Motorized Lateral Drive
- Lateral Travel: (± 4 ") ± 100 mm

Accessories



Wireless Hand Control



Foot Control for all Movements



Accessory Clamp



Patient Arm Board and Pad

Specifications

0047-ELITE SERIES SPECIFICATIONS		TM-5	TM-4	TM-3	TM-2	TL-2	TT-2	TE-1
(1) Motorized Elevation (Lift) Travel	Min.	31" 790mm	30" 765mm	29" 745mm	23" 580mm	26" 650mm	26" 650mm	22" 560mm
300mm Lift Standard	Max.	43" 1090mm	42" 1065mm	41" 1045mm	35" 890mm	37" 950mm	37" 950mm	34" 860mm
(1) Motorized Elevation (Lift) Travel	Min.	NO	NO	NO	25" 635mm	28" 710mm	28" 710mm	24" 610mm
400mm Lift (Optional)	Max.	NO	NO	NO	41" 1040mm	43" 1110mm	43" 1110mm	40" 1010mm
(2) Motorized Trendelenburg Tilt		± 15°	± 15°	± 15°	NO	NO	± 15°	NO
(3) Motorized Longitudinal Tabletop Travel		15.75" 400mm	15.75" 400mm	NO	15.75" 400mm	NO	NO	NO
(4) Motorized Lateral Tabletop Tilt		± 15°	± 15°	± 15°	NO	± 15°	NO	NO
(5) Motorized Lateral Tabletop Travel		± 4" 100mm	NO	NO	NO	NO	NO	NO
Patient Weight Capacity		500lb 227kg	500lb 227kg	500lb 227kg	500lb 227kg	500lb 227kg	500lb 227kg	500lb 227kg
Carbon Fiber Tabletop 22" wide x 82" long, with Facial access cutout		YES	YES	YES	YES	YES	YES	YES
Carbon Fiber Rectanaglar Tabletop 22" wide x 82" long, with Facial access cutout		OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL
Unobstructed Tabletop Imaging Area		58" 1473mm	58" 1473mm	58" 1473mm	58" 1473mm	58" 1473mm	58" 1473mm	58" 1473mm
Radiolucent Sofcare Pad		YES	YES	YES	YES	YES	YES	YES
Pateint Restraints (3)		YES	YES	YES	YES	YES	YES	YES
OPTIONS								
Dual Patient Arm Boards with Pads		OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL
Snap-on Arm Board with Pad and Mounting Clamp		OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL
Wireless Hand Control with Receiver		NO	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL
Foot Pedal Control		OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL
Rectangular Tabletop		OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL



Technical specifications are subject to change without notice.
 © 2009 Arcoma North America, Inc 0047-005-A (03-09)



North American Office 23151 Alcalde Dr., Suite C-8, Laguna Hills, CA 92653
 Worldwide Headquarters Anavagen 1, Vaxjo SE-325 46, Sweden